PHILIPS

SCIENTIFIC EQUIPMENT

AREOMETER PW 7020

DIRECTIONS FOR USE

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MODE D'EMPLOI

Areometer for determining the purity of liquid nitrogen contaminated with oxygen (type PW 7020)

The density of liquid nitrogen depends on the oxygen content and the boiling point, while the latter quantity is affected by the barometric pressure.

From the above it appears that — after a reading of the barometer has been taken — the purity of liquid nitrogen may be determined by measuring the density with the aid of an areometer.

Admixtures of argon also affect the degree of purity, but within the measuring range of 99-100 % N₂ the liquid will rarely contain argon. Nevertheless, should the liquid contain any argon, a greater degree of contamination is measured than corresponds with the facts.

The measuring accuracy is better than 0.2 %.

Always use the graph of the same number as that of the areometer

1. Clean the areometer and place it into the measuring vessel.
2. Fill the vessel carefully with the liquid nitrogen by pouring it from a low level.
3. Read the areometer after the liquid has come to rest.
4. Take down the barometric pressure.
5. Read the purity on the relevant graph.
6. Replace everything in the case after use; this prevents damage.

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